

Product Safety Information

Safety cabinet type 90 asecos Q-CLASSIC-90, hinged doors, W 1193 x D 615 x H 1953 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 461739, 461738, 461737, 461736, 461735, 461734, 461733

Safety Instructions:

- Only use the safety storage cabinet to store toxic, flammable and water-polluting substances in accordance with the regulations.
- Always secure the cabinet with the profile cylinder lock to prevent unauthorised access. Keep the key in a safe place.
- Ensure that the cabinet has properly ventilated exhaust air connections to remove hazardous vapours and ensure an explosion-proof environment.
- Avoid exceeding the maximum load of 50 kg for the tray shelves.
- Keep the cabinet in a dry and well-ventilated area to minimise corrosion and the risk of vapours or leaks.
- Observe the local regulations for the storage of hazardous substances.
- Make sure that the cabinet is stable and cannot tip over, especially when loading and unloading.
- The safety storage cabinet should be checked regularly for functionality and safety.
- Provide technical documentation and safety information in an understandable language. Enclose proof of conformity with European standards.
- Only use the cabinet indoors and connect it to a technical exhaust air system.
- Fire resistant (F90): Avoid open flames and excessive heat in the vicinity of the cabinet.
- Ensure regular maintenance of the ventilation systems to ensure safe operation.
- Use floor drip pans to prevent leaks and contamination.
- Read the operating instructions thoroughly before use and ensure that the storage areas are safe and well ventilated.
- If there is visible damage, do not use the cabinet and contact the manufacturer.