

Product Safety Information

Laboratory balance KERN EMB

Manufacturer/Responsible Person: Kern & Sohn GmbH, Ziegelei 1, 72336 Balingen, Deutschland

E-Mail: info@kern-sohn.com

Please note the technical data of the products!

Article Numbers: 67244, 67243, 157215, 157214

Safety Instructions:

- Batteries must be kept out of the reach of children.
- Do not use on an uneven surface or unstable surfaces. Use on a stable, level surface.
- Avoid pouring liquids onto the appliance or using it in a damp environment - risk of short circuit! Do not immerse in water to avoid possible electrical damage.
- Read the operating instructions carefully before use.
- Do not place the appliance near heat sources.
- Use only for the purposes specified in the manual. Use only within the intended weighing range: Excess weight can cause damage. Do not overload: Do not exceed the maximum weight of 5200 g.
- Small parts (even in the packaging) can be swallowed - choking hazard! Risk of injury to children, therefore keep out of reach. Not suitable for children under 3 years of age without supervision.
- Observe safety precautions: Avoid risk of falling or breakage, ensure safe handling.
- Do not overheat batteries or throw them into a fire - fire hazard! Dispose of batteries in accordance with the manufacturer's instructions.
- The batteries may only be changed by adults.
- Keep hands dry: Wet hands can impair consistency and readability.