

Product Safety Information

Lab swivel chair Labster 9103

Please note the technical data of the products!

Article Numbers: 141364, 141362, 141360

Safety Instructions:

- Use only in suitable laboratory environments to minimise safety risks. The chair must not be used in public areas as it is intended for specialised medical activities.
- Regularly check for mechanical damage, especially the seat and backrest adjustment, as well as the mechanical functions of the chair to ensure safe working.
- Ensure that the chair is used on stable floors to prevent it from tipping over.
- Avoid excessive load and weight on the chair to ensure the durability and ergonomic properties of the components.
- Checking and regular maintenance of materials with suitable cleaning agents and disinfectants to maintain hygiene in laboratory environments.
- Read the instructions for use and safety information carefully before first use.
- When using the laboratory swivel chair, ensure the correct height adjustment to prevent falls. The height adjustment is designed between 450 - 650 mm; use outside this range may result in injury. Secure the height adjustment properly to prevent it from slipping.
- It is recommended that only authorised disinfectants are used for cleaning in order to prevent damage to the material.
- The chair is designed so that it can be used in hygiene-critical areas. Hygiene requirements should be observed.
- Please ensure that the chair is used correctly in accordance with the manufacturer's instructions.
- Avoid using the product in areas with extreme temperatures or high humidity.