

Product Safety Information

Standing aids Futura Series AF6

Manufacturer/Responsible Person: Mey Systems GmbH, Merlach 16, 96145 Sesslach-Merlach, Deutschland

E-Mail: info@mey-components.de

Please note the technical data of the products!

Article Numbers: 45571, 45570, 45569

Safety Instructions:

- Use only in accordance with the operating instructions.
- Suitable for people with a maximum weight of 120 kg.
- Do not place excessive weight on the seat to prevent damage or accidents.
- Check the product regularly for damage or signs of wear, especially before it is used.
- Ensure that the seat mechanism works properly to prevent falls.
- Adjust the height setting correctly and safely and ensure that the height adjustment is locked to prevent slipping during use.
- Observe safety precautions when adjusting the height.
- Check that the standing aid is used on a stable and level surface to prevent it from tipping over.
- Users should be advised not to use the standing aid if they have physical limitations or are unsteady in their balance.