

Product Safety Information

Chemical cabinet Asecos CS-CLASSIC, hinged door, 2 pull-out trays, W 545 x D 520 x H 1105 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 132155, 132154, 132153, 132152, 132151, 132150, 132149

Safety Instructions:

- Use the product only for the storage of water-polluting, non-flammable liquids and toxic substances in accordance with regulations and not for flammable materials.
- Ensure that the cabinet is installed in a well-ventilated area, is naturally and technically ventilated and that the ventilation grilles are not blocked.
- The cabinet should only be used by trained personnel for the intended purpose in order to minimise the risk of spillage of hazardous substances.
- Ensure that the cabinet is securely locked and that the cylinder lock is used to prevent access by unauthorised persons.
- Do not allow children access to the stored materials and avoid potential choking hazards due to small parts.
- Check the cabinet regularly for damage or wear and tear and check the plastic drip trays for damage.
- Keep the cabinet away from direct sunlight, moisture and at a safe distance from heat sources.
- Provide safety instructions and instructions in an understandable language and comply with local storage requirements.
- In the event of damage or faults, contact the manufacturer immediately.
- Observe environmental regulations when disposing of chemicals.
- Ensure secure installation of the cabinet, for example by using the levelling feet for stable positioning.