

Product Safety Information

Chemical cabinet Asecos CS-CLASSIC-G, hinged door and glass cut-out, 2 shelves, H 1105 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 132148, 132147, 132146, 132145, 132144, 132143, 132142

Safety Instructions:

- Take precautionary measures when handling water-polluting substances; use a cabinet to store water-polluting, non-flammable liquids and toxic substances in accordance with regulations.
- Store non-flammable liquids safely and observe precautionary measures; suitable for water-polluting liquids, but not for flammable substances.
- Ensure regular inspection and proper installation of ventilation systems; natural ventilation grilles must not be blocked, keep ventilation grilles near the floor free.
- Protect the cabinet from improper use and ensure stable installation; always keep the doors closed and use a cylinder lock for security.
- Follow protective measures and wear protective equipment when handling toxic substances and chemicals.
- Observe restrictions for use by small children; place the cabinet out of the reach of children.
- In the event of an accident, notify the relevant authorities immediately and contact the manufacturer in the event of damage or leaks.
- Stably adjust height-adjustable shelves to specification to ensure stability and utilise flexible storage options.
- Avoid overfilling the cupboard to ensure safer handling.
- Check the cabinet regularly for damage and leak-tightness of the base drip tray, replace if necessary.
- Always keep the door closed during operation to minimise the risk of vapours or outgassing.
- Take protective measures when handling and storing toxic substances, users should wear appropriate protective equipment.
- Recall measures must be observed immediately in the event of product safety risks.
- Before connecting to technical ventilation, ensure that it functions properly; in the case of flammable liquids, the cabinet must be connected to such ventilation.
- Observe and keep all safety instructions and product information in understandable language; failure to do so may result in leaks.
- Not suitable for the storage of flammable substances; ensure use only in well-ventilated areas.
- Suitable for the storage of water-polluting, non-flammable liquids and toxic substances in accordance with regulations; exercise caution during storage.