

Product Safety Information

Chemical cabinet Asecos CS-CLASSIC-G, hinged door w. glass cut-out, 3 shelves, W 1055 x D 520 x H 1950 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 132064, 132063, 132062, 132061, 132060, 132059, 132058

Safety Instructions:

- Use for the storage of water-polluting and non-flammable liquids and toxic substances in work areas in accordance with regulations. The cabinet is not suitable for flammable liquids.
- Access to the cabinet is only permitted to authorised personnel. Children and unauthorised persons must be kept away.
- Ensure that the cabinet is always locked when not in use; this is done using the integrated cylinder lock.
- Please observe the ventilation requirements. The cabinet is prepared for connection to a technical ventilation system and regular checks are necessary to prevent the accumulation of hazardous vapours.
- Safety information and instructions for use must be clearly visible and provided in understandable language, in accordance with applicable standards.
- Ensure that the cabinet is stabilised against tipping over and is installed on a level surface. Carry out regular checks for signs of corrosion or damage.
- Empty containers must be properly labelled and disposed of safely.
- Use the cabinet in accordance with the manufacturer's guidelines and only for its intended purpose to minimise potential risks. Do not mix with other chemicals.
- If the product or the locking mechanisms are damaged, take immediate action to repair or replace them.
- Protect the cabinet from environmental influences and direct sunlight.
- Use ventilation grilles near the floor for natural ventilation.
- Observe the operating instructions when using and installing the cabinet.