

Product Safety Information

Chemical cabinet Asecos CS-CLASSIC-G, hinged door w. glass cut-out, 2 pull-out trays, W 1055 x D 520 x H 1105 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 132113, 132112, 132111, 132110, 132109, 132108, 132107

Safety Instructions:

- Safe storage of water-polluting, non-flammable liquids and toxic substances in accordance with the applicable regulations. Not suitable for flammable liquids.
- The cabinet must be checked regularly for damage, leaks and rust.
- Before using the appliance, ensure that it is installed correctly, that all safety functions are working and that you have read the user manual.
- Ventilation must not be obstructed. Connect technical ventilation if necessary to avoid overpressure or hazardous vapours.
- Do not overload the cabinet or use it improperly to avoid overloading and possible accidents.
- Lock the cabinet securely with a cylinder lock to prevent unauthorised access. Keep doors closed.
- Only use the product in well-ventilated areas. Do not block ventilation grilles to ensure natural ventilation.
- There is a risk of chemical accidents if used improperly. Follow instructions for the safe use and disposal of hazardous materials.
- User instructions and safety instructions must be provided clearly and comprehensibly. For commercial use only, not suitable for private households.
- Ensure that the drip trays and height-adjustable shelves are intact and securely fastened.
- Choose a suitable location that is protected from direct sunlight and moisture.
- In the event of faults or damage to the cabinet, contact the manufacturer immediately.
- For added safety, the product should be coloured safety yellow and protective equipment should be worn when handling potentially hazardous substances.