

Product Safety Information

Chemical cabinet Asecos CS-CLASSIC-G, hinged door w. glass cut-out, 4 pull-out trays, W 1055 x D 520 x H 1950 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 132057, 132056, 132055, 132054, 132053, 132052, 132051

Safety Instructions:

- Only store highly hazardous liquids in well ventilated areas and ensure that the cabinet is properly ventilated to prevent dangerous vapour formation.
- Ensure that safety markings and instructions are always followed.
- Store and use cabinets for the storage of toxic substances, water-polluting and non-flammable liquids, and avoid using them for flammable liquids.
- Ensure that the cabinet is kept out of the reach of children and is protected with a functioning cylinder lock to prevent unauthorised access.
- Regularly check the lock, the tightness, stability and ventilation of the structure.
- Use personal protective equipment (e.g. gloves, safety goggles) when handling hazardous chemicals.
- When arranging chemicals in the cabinet, they should be arranged according to compatibility in order to prevent dangerous reactions.
- Keep a certain distance from heat sources and protect the cabinet from direct sunlight to prevent material fatigue.
- Ensure that the cabinet is installed on a firm and level floor to ensure stability and avoid blocking ventilation grilles.
- Cabinets should only be used by specialised personnel; sale only to specialised personnel or companies that handle hazardous substances.
- Ensure that the cabinet is correctly assembled and secured and that the shelves and drip trays are properly installed to prevent leakage.
- Avoid combining with other chemical products to prevent dangerous chemical reactions.
- Improper use may result in the release of toxic vapours or liquids.