

Product Safety Information

LABSIT laboratory swivel chair, without armrest, with castors, backrest height 420 mm

Please note the technical data of the products!

Article Numbers: 144311, 144310, 144309, 144308, 144307, 144306, 144305, 144304, 144303, 144302

Safety Instructions:

- Make sure that the chair is used on a stable, hard surface to prevent it from tipping over.
- The seat height should only be continuously adjusted using the safety gas spring in accordance with the instructions in order to avoid the risk of injury due to malfunctions.
- Before use, check the chair for visible damage, especially to the gas spring, and do not use it if there are any defects.
- Only use the castors on suitable hard floors; avoid uneven or soft floors to prevent stability problems.
- Follow the enclosed assembly instructions carefully to ensure that the chair is assembled correctly and falls are avoided.
- Keep the castors and seat clean for safe working and avoid use on slippery floors without suitable mats.
- The chair is not suitable for use by children and should not be used outdoors.
- Regular checks of the stability and functionality of all parts, especially the gas spring and castors, are essential to ensure safe use.
- Do not exceed the maximum load capacity of the chair to avoid damage.
- Ensure correct maintenance of the washable artificial leather and fabric covers to guarantee hygiene and safety.