

Product Safety Information

Safety storage cabinet type 30 asecos Q-CLASSIC-30, hinged door, stop left/right, W 564 x D 620 x H 1947 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 462229, 462228, 462227, 462226, 462225, 462224, 462223, 462222, 462221, 462220, 462219, 462218, 462217, 462216

Safety Instructions:

- Ensure that the safety storage cabinet is properly installed, anchored and stable to prevent it from tipping over.
- Ensure compliance with the applicable regulations and safety data sheets for the storage of hazardous chemicals and substances.
- Use personal protective equipment when handling stored substances.
- Do not leave the door of the cupboard unsecured and always lock it with a profile cylinder lock to prevent unauthorised access.
- The cabinet must be installed in a well-ventilated place and connected to a technical exhaust air system to avoid vapours and possible explosion hazards.
- Carry out regular maintenance and inspection of the ventilation system and the tightness of the drip tray.
- Do not exceed the maximum load capacity of the shelves: observe the load limit of 75 kg per shelf.
- Only use for the safe storage of GHS categories 1-3 flammable, toxic or aggressive liquids.
- The cabinet should not be located near sources of ignition, heat sources or direct sunlight.
- Fire-retardant (F30), but not fireproof: Protect from direct ignition and high temperatures.
- Store chemicals and substances only for their intended use and in accordance with the safety data sheets.
- Clean hands thoroughly after use to avoid chemical residues.
- Regular inspection of the integrated ventilation system required, ensure that exhaust air connections are not blocked.
- Only allow authorised persons access to the cabinet.
- Read and follow the instructions and safety information carefully before use.