

Product Safety Information

Safety cabinet type 90 asecos S-CLASSIC-90, hinged door, stop left/right, W 596 x D 616 x H 1968 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 462334, 462333, 462332, 462331

Safety Instructions:

- Observe fire resistance F90 and use a fire-resistant cabinet in accordance with the F90 standard.
- Safety yellow front for better visibility of chemicals.
- Always keep doors closed to minimise chemical vapours and the spread of fire and to prevent unauthorised access.
- Only use the door hold-open system under safe conditions.
- Protect contents from unauthorised access - insert profile cylinder lock and check regularly that it is working properly.
- Use a technical exhaust air system for explosion protection and ensure that the enclosure is connected to a wall ventilation system.
- Ensure an air exchange rate of 10 times and carry out regular maintenance and inspection of the ventilation system.
- Store chemicals in accordance with GHS regulations and store safely in accordance with regulations for toxic and aggressive chemicals.
- Wear suitable personal protective equipment when working with toxic or hazardous substances.
- Only use indoors and ensure that the cabinet is not used to store unauthorised substances.
- Do not use naked flames or ignition sources near the cabinet.
- Ensure the correct disposal of hazardous waste from the cabinet.
- Never overfill the cupboard (observe the maximum load capacity).
- Attention: Only suitable for the storage of toxic and flammable substances.
- Observe safety precautions: Always close the door lock securely.