

Product Safety Information

Safety cabinet type 90 asecos S-CLASSIC-90, hinged door with door locking system, left/right hinges, W 596 x D 616 x H 1968 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 462350, 462349, 462348, 462347

Safety Instructions:

- Store flammable liquids safely in accordance with applicable regulations for toxic and aggressive chemicals.
- For indoor use only.
- Always lock the door with the profile cylinder lock and keep the cabinet locked to prevent unauthorised access.
- A proper exhaust air system is required for operation to avoid the risk of explosion. Check regularly and ensure that the ventilation system is not blocked.
- Observe explosion protection when storing hazardous substances.
- Carry out regular visual inspections for possible damage and check the leak tightness of the drip trays to prevent leaks.
- Wear suitable protective equipment when operating the cabinet and avoid direct contact with chemical substances.
- Only use compatible chemicals that are suitable for storage in this cabinet and store the chemicals authorised for the cabinet, observe maximum filling quantities.
- Read the operating instructions carefully before use.
- Use the door hold-open device to prevent the door from closing unintentionally and to avoid the risk of jamming when closing the door.
- The product is fire-resistant, but not protected against open flames.
- Follow the applicable regulations for the disposal of hazardous waste and chemicals.
- Ensure that the cabinet is installed in a well-ventilated area to avoid hazardous vapours.
- In the event of chemical accidents, immediately take appropriate measures to avert danger and inform the responsible authorities.
- Only use the safety storage cabinet in accordance with the manufacturer's instructions.
- Do not exceed the load capacity of the pull-out trays and floor drip tray.