

Product Safety Information

Safety base cabinet type 90 asecos UB-T-90, hinged doors, W 1102 x D 502 x H 600 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 462833, 462832

Safety Instructions:

- Safe storage of toxic or aggressive chemicals and flammable liquids in accordance with the GHS guidelines.
- Only authorised persons are permitted to use the cabinet and ensure that it is only used for approved liquids and chemicals.
- Regularly check the ventilation system for proper function; use it in accordance with the instructions to avoid explosive atmospheres.
- Check cylinder lock regularly for accessibility and security, keep key out of reach of unauthorised persons.
- Do not block the hinged doors to ensure safe access.
- The cabinet must not be overloaded, observe the maximum load capacity. Take into account the load capacity of the perforated plate insert (75 kg) and the floor drip tray (22 litres).
- Floor drip tray can only be used with drained chemicals.
- Safe storage of safety data sheets (SDS) for all stored chemicals in order to be able to act quickly in an emergency.