

Product Safety Information

Chemical cabinet asecos CS-CLASSIC, hinged doors, W 810 x D 520 x H 1950 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 463037, 463036, 463035, 463034, 463033, 463032, 463031, 463030, 463029, 463028, 463027, 463026, 463025, 463024

Safety Instructions:

- For indoor use only.
- Do not store with flammable or explosive substances; use special chemical cabinets.
- Provide safety instructions in understandable language and operating instructions clearly.
- Keep out of reach of children.
- Ensure regular maintenance and inspection of the ventilation system.
- Observe the maximum load capacity of the shelves (60 kg).
- Store chemicals only in suitable containers recommended by the manufacturer.
- Use hinged doors and a cylinder lock to prevent unauthorised access.
- Connect the cabinet to a technical exhaust air system with an exhaust air connection (75 mm) and have it installed by qualified personnel only.
- Use of personal protective equipment required when handling hazardous substances.
- Check the floor drip tray regularly for leaks and condition.
- Ensure safe storage of hazardous substances in accordance with legal regulations.
- Only use the chemical cabinet for its intended purpose and check it regularly for structural integrity and functionality.
- Ensure that the chemical cabinet is explosion-proof and properly ventilated.
- Dispose of chemicals safely in accordance with local regulations.
- Secure the cabinet against tipping over and protect it from direct sunlight and extreme temperatures.
- In the event of chemical releases or leaks, have emergency measures ready.