

Product Safety Information

Chemical cabinet asecos CS-CLASSIC, hinged doors, W 810 x D 520 x H 1105 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 463084, 463083, 463082, 463081, 463080, 463079, 463078, 463077, 463076, 463075, 463074, 463073, 463072

Safety Instructions:

- Use only for the storage of hazardous substances and ensure that they are stored in accordance with legal regulations.
- Avoid overloading the shelves (max. load capacity 60 kg) and only use the floor drip tray for its intended purpose (max. 33 litres).
- For indoor use only; protect the cabinet from moisture, direct sunlight and excessive heating.
- Ensure that the cabinet is always well ventilated and properly operated with a technical exhaust air system.
- Carefully lock hinged doors with a cylinder lock to prevent unauthorised access.
- Access only for authorised persons and only to be operated by trained personnel.
- Regularly check the cabinet and its components (doors, locks, ventilation system) for damage and inform the manufacturer if any damage is found.
- Use personal protective equipment (PPE) such as gloves and safety goggles when handling hazardous substances.
- Keep the hinged doors closed when the cabinet is in use and lock it when not in use.
- In the event of accidents or leaks, take immediate action, close the cabinet and inform the relevant authorities.
- Install a suitable exhaust air system for explosion protection and ensure regular inspection of the ventilation system and exhaust systems.
- Ensure safe storage of chemical substances in compliance with safety regulations and safety data sheets.
- Ensure compliance with the safety regulations for hazardous substances in operation and access to an extinguishing agent in the vicinity.