

Product Safety Information

Chemical cabinet asecos CS-CLASSIC, hinged door, stop left, W 545 x D 520 x H 1950 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 463141, 463140, 463139, 463138, 463137, 463136, 463135, 463134, 463133, 463132, 463131, 463130, 463129, 463128

Safety Instructions:

- Read and follow the safety instructions and operating instructions carefully.
- Store and handle chemicals in accordance with local regulations and safety data sheets.
- For indoor use only. The chemical cabinet is not suitable for outdoor use.
- Ensure that the cabinet is properly ventilated and check the ventilation system regularly.
- Max. Load capacity of the shelves must not exceed 50 kg. Keep the cabinet in an upright position.
- Always lock the cabinet when not in use to prevent unauthorised access. Use the cylinder lock.
- Only authorised persons have access to the cabinet. Keep out of the reach of children.
- The cabinet may only be used for the storage of authorised hazardous substances. Security against unauthorised access must be ensured.
- Explosion hazard: Keep the cabinet away from sources of ignition. Use a technical exhaust air system for explosive chemicals.
- In the event of damage or leaks, react immediately and inform the manufacturer or supplier.
- Dispose of packaging and chemicals safely in accordance with the regulations. Observe the return regulations.
- Monitor the ventilation and exhaust air connection (recommended: 75 mm) regularly to minimise risks.
- Handle door with glass cut-out carefully to avoid breakage.
- Wear PPE when handling hazardous substances and alert the fire brigade and evacuate immediately in the event of a fire.
- Regularly check the cabinet for damage, locking mechanisms and exhaust air systems.