

Product Safety Information

Chemical cabinet asecos CS-CLASSIC, hinged door, stop left, W 545 x D 520 x H 1105 mm

Manufacturer/Responsible Person: asecos GmbH, Weiherfeldsiedlung 16-18, 63584 Gründau, Deutschland

E-Mail: info@asecos.com

Please note the technical data of the products!

Article Numbers: 463160, 463159, 463158, 463157, 463156, 463155, 463154, 463153, 463152, 463151, 463150, 463149

Safety Instructions:

- Use in accordance with the requirements for the safe handling of hazardous substances.
- Special care must be taken when opening the hinged door to prevent the release of hazardous substances.
- Maximum load capacity of the shelves: 50 kg. Do not overload to avoid the risk of falling.
- The exhaust air connection must be properly connected and maintained to avoid the risk of explosion.
- Place the chemical cabinet out of the reach of children and unauthorised persons.
- Wear suitable personal protective equipment (PPE) when working with chemicals.
- Only use the product indoors and protect it from direct sunlight and extreme heat.
- Observe explosion protection: Connection to a technical exhaust air system required.
- Do not store any chemical substances outside the designated areas.
- Carry out regular inspection and maintenance of the oven and the ventilation system.
- Handle hinged doors with care and secure with a cylinder lock to prevent injury and unauthorised access.
- Use a suitable cabinet for non-flammable chemicals. Do not store flammable chemicals.
- Keep safety data sheets available for all stored chemicals and update them regularly.
- In the event of leaks: Use a floor drip tray with a volume of at least 22 litres and check regularly.
- Ensure that all employees are appropriately trained.