

Product Safety Information

Safety check for HeartSine defibrillators, according to MPBetreibV § 11, incl. loaner device

Manufacturer/Responsible Person: MedX5 GmbH & Co. KG, Humboldtstr. 9, 86316 Friedberg, Deutschland

E-Mail: info@medx5.com

Please note the technical data of the products!

Article Number: 178314

Safety Instructions:

- Ensure that defibrillators are always in working order
- Avoidance of electric shocks due to improper use
- Carry out regular checks in accordance with the applicable ERC & AHA guidelines
- Document test reports and protocols
- Use of test devices authorised by the manufacturer
- Replacing batteries and updating software as required
- Clean and disinfect appliances regularly
- Ensure the use of replacement devices during maintenance
- Visual inspections must not replace technical inspections
- Do not accept any unofficial test protocols