

Product Safety Information

Standing aid meychair Futura Light AF-SR, folding & seat height adjustable, without step-up aid, steel, black, various variants upholstery

Manufacturer/Responsible Person: Mey Systems GmbH, Merlach 16, 96145 Sesslach-Merlach, Deutschland

E-Mail: info@mey-components.de

Please note the technical data of the products!

Article Numbers: 143022, 143020, 143018

Safety Instructions:

- Do not exceed the maximum load capacity of 120 kg.
- Observe the safety instructions for installation and use.
- Before use, ensure that all mechanical parts are functioning correctly.
- Do not use as a climbing aid as it is not designed for this purpose.
- If used incorrectly, it may tip over or cause injury.
- Place on a stable surface before use to prevent tipping or falling over.
- Regularly check for damage or wear to the components, especially the height adjustment.
- Do not adjust the seat height abruptly to avoid injury from the moving parts.
- Always use the product in accordance with the manufacturer's instructions.
- Check the seat height to ensure that you are sitting safely.
- Use the space-saving folding function to avoid injuries.