

Product Safety Information

Ergotron CareFit™ 2.0 documentation trolley, with hinge for display, up to 27", 2 storage levels, on castors, white/warm gray

Manufacturer/Responsible Person: ALSO Deutschland GmbH, Lange Wende 43, 59494 Soest, Deutschland

E-Mail: info@also.com

Please note the technical data of the products!

Article Number: 660894

Safety Instructions:

- Risk of suffocation due to improper use of the storage levels.
- Risk of tipping over if the trolley is used incorrectly or overloaded.
- When using electrical appliances on the trolley, cable management must be ensured to avoid the risk of tripping.
- The trolley should not be used in damp areas in order to minimise electrical risks.
- Unplug appliances before folding or moving the trolley.
- Observe the maximum load capacity of 9.1 kg to avoid damage or injury.