

Product Safety Information

Laboratory balance Kern EMS 6K0.1/12K0.1, readability 0.1 g, weighing range 6/12 kg, large LC display, W 200 x D 280 x H 65 mm, plastic, white

Manufacturer/Responsible Person: Kern & Sohn GmbH, Ziegelei 1, 72336 Balingen, Deutschland
E-Mail: info@kern-sohn.com

Please note the technical data of the products!

Article Numbers: 131409, 131408

Safety Instructions:

- Ensure that the scale is used within a temperature range of +5 °C to +35 °C.
- When using in battery mode, ensure that the batteries are inserted correctly.
- The mains power supply must not be overloaded in order to avoid the risk of fire.
- The laboratory balance must not be used in a damp environment in order to avoid electrical damage.
- The smallest parts or the scales themselves must not be used or reached by children in order to avoid the risk of suffocation.
- Clean the surfaces of the scales regularly so as not to impair the functions.
- For indoor use only.
- Electronic components are not waterproof. Protect from moisture.
- Do not overload. Observe the maximum weighing range of 12 kg.
- Do not use on uneven surfaces to avoid incorrect weights.
- In the event of faults or damage, disconnect from the mains immediately and contact customer service.