

Product Safety Information

Laboratory balance Kern EMS 6K1/12K1, readability 1 g, weighing range 6/12 kg, large LC display, W 200 x D 280 x H 65 mm, plastic, white

Manufacturer/Responsible Person: Kern & Sohn GmbH, Ziegelei 1, 72336 Balingen, Deutschland
E-Mail: info@kern-sohn.com

Please note the technical data of the products!

Article Numbers: 131467, 131407

Safety Instructions:

- The laboratory balance must not be used in damp or wet environments.
- Make sure that the device is placed on a stable and level surface to ensure accurate measurement and to prevent it from tipping over.
- The laboratory balance is not suitable for use by children.
- Avoid overloading the appliance; the maximum capacity is 12 kg.
- Protect the laboratory balance from extreme temperatures above +35 °C or below +5 °C.
- Avoid attaching external devices that could influence the function of the scales.
- Check the battery status regularly; use batteries at your own risk to avoid malfunctions.
- Do not use the laboratory balance improperly to avoid damage to electrical components.
- To avoid contamination, do not use the laboratory balance in the vicinity of microbes or chemicals.
- Calibrate the laboratory balance regularly to ensure precise measurements.