

Product Safety Information

KERN transmitted light microscope OBE 121

Manufacturer/Responsible Person: Kern & Sohn GmbH, Ziegelei 1, 72336 Balingen, Deutschland

E-Mail: info@kern-sohn.com

Please note the technical data of the products!

Article Numbers: 662481, 662480, 662479

Safety Instructions:

- Ensure that the microscope is stable during use to prevent falls and damage.
- Keep the microscope and its components away from water and damp environments to avoid electrical hazards.
- Store the microscope in a dry place to prevent damage from dust or moisture.
- Avoid direct contact with the LEDs during operation as they can become hot and there is a risk of burns.
- Children should only use the microscope under supervision to avoid injury and the risk of choking on small accessories.
- Ensure that the power supply is connected correctly to avoid overheating and electric shocks.
- Avoid direct light on the eyes to prevent eye damage from the LED.
- Dispose of electronic components of the microscope in an environmentally friendly manner in accordance with local regulations.
- Only use accessories and components recommended by the manufacturer.