

Product Safety Information

Laboratory swivel chair LEANOS 9133, imitation leather, Stamskin Top blue

Please note the technical data of the products!

Article Numbers: 41267, 41266, 41265, 307992, 307991, 307990, 307989

Safety Instructions:

- Ensure that the maximum weight load of the chair is not exceeded in order to avoid mechanical damage and overloading.
- Avoid excessive tilting of the chair and ensure that it is used on a stable and level surface to prevent tipping.
- Correctly adjust the seat height and tilt before use to avoid injury.
- Regularly check the stability and functionality of the castors and the gas spring to minimise the risk of falling.
- Do not use aggressive cleaning agents when cleaning to avoid damaging the material.
- Safety gas spring must be correctly positioned and maintained.
- Beware of choking hazard due to small parts that can be swallowed or loose parts.
- Be aware of the risk of slipping due to improper use, especially on slippery or moving surfaces.
- Product not suitable for use in high-risk areas or for children and vulnerable persons.
- Avoid using the chair if it is damaged.
- Regularly check the stability of the chair and keep fingers away from moving parts.