

Product Safety Information

Laboratory swivel chair LEANOS 9138, imitation leather, Stamskin Top blue

Please note the technical data of the products!

Article Numbers: 41292, 41291, 41290, 308025, 308024, 308023, 308022, 308021, 308020

Safety Instructions:

- Ensure that the chair height is set correctly to avoid the risk of falling and injury due to improper use of the gas spring.
- Keep the chair clean and clean it regularly to prevent bacteria build-up and material damage.
- Check the product regularly for wear, damage and proper functioning of the castors to ensure safety.
- Observe and comply with the manufacturer's instructions and safety information to ensure safe use.
- Avoid risk of suffocation from loose or small moving parts, use only stabilised.
- Avoid the risk of slipping and tipping when using on uneven or unsuitable floors, only use on level and stable floors.
- Intended for adult use only to minimise the risk of injury.
- Do not overload the chair to prevent it from breaking or tipping over; observe the maximum utilisation load of 120 kg.
- Avoid moisture and extreme temperatures to protect the material.
- Ensure correct installation and maintenance of the safety gas spring in accordance with the manufacturer's instructions.
- Regularly check the gas spring for safety, replace if damaged.
- Provide the chair with clear and understandable safety instructions and warn users of hazards due to improper use.
- Only use in recommended environments such as laboratories or workshops.
- Use suitable cleaning agents for cleaning and care to prevent damage to the artificial leather.